

PE Anti-Human CD133(W6B3C1)

Product Information

Catalog No.:	PE-30154
Conjugation:	PE
Size:	20T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	W6B3C1
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	Prominin-1 AC133 Antigen AC133 CD133 CORD CORD12 Hematopoietic stem cell antigen MCDR MCDR2 MSTP MSTP061 PROM PROM1 PROML PROML1 Pr ominin-1 RP41 STGD STGD4 hProminin
Uniprot ID :	O43490

Application

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image: