

PE Anti-Human CD66b (G10F5)

Product Information

Catalog No.:	PE-30156
Conjugation:	PE
Size:	20T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	G10F5
IsoType:	IgGM, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	CD66b CD67 CD67 Antigen CEACAM CEACAM8 CGM CGM6 Carcinoembryonic Antigen CGM Carcinoembryonic Antigen CGM6 Carcinoembryonic antigen-related cell adhesion molecule 8 NCA NCA-95 Non-Specific Cross-Reacting Antigen NCA
Uniprot ID :	P31997

Application

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image: