

PE-FineTest®594 Anti-Mouse CD4 (GK1.5)

Product Information

Catalog No.:	PE5-30001
Conjugation:	PE-FineTest®594
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	GK1.5
IsoType:	IgG2b, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	T-cell surface glycoprotein CD4 T-cell surface antigen T4 Leu-3 L3T4 T4
Uniprot ID :	P06332

Application

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



C57BL/6 murine splenocytes are stained with PE/FineTest® 594 Anti-Mouse CD4 Antibody (filled red histogram). Unstained splenocytes (empty black histogram) are used as control.