
Biotin Anti-CPS1

Product Information

Catalog No.:	biotin-01934
Conjugation:	biotin
Size:	100µg
Form:	liquid
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 5mg/ml BSA, 50% glycerol, pH7.6. -20°C for 24 months, avoid repeated freezing and thawing.

Immunogen information

Synonyms:	Carbamoyl-phosphate synthase [ammonia], mitochondrial Carbamoyl-phosphate synthetase I (CPSase I) CPS1
Uniprot ID :	P31327

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA