

## FITC conjugated Affinipure protein AG

<b>Catalogue No.:</b>	FNSA-0132
<b>Species:</b>	Staphylococcus aureus
<b>Uniprot ID:</b>	P38507
<b>Host:</b>	E.Coli
<b>Tags:</b>	N-terminal His Tag
<b>Purity:</b>	Greater than 90% as determined by SDS-PAGE
<b>Formulation:</b>	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, pH7.6
<b>Storage:</b>	2-8°C for 12 months
<b>Synonyms:</b>	proteinAG
<b>Safety Note:</b>	<p>This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.</p>