


Recombinant Human GABRG2-2

Catalogue No.:	P2078
Species:	Human
Uniprot ID:	P18507
Expression Region:	361-443
Host:	E.Coli
Tags:	N-terminal His-IF2DI Tag
Molecular Weight:	30.3 kDa under reducing conditions
Purity:	Greater than 90% as determined by SDS-PAGE.
Formulation:	Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 150 mM NaCl with 5% trehalose, pH7.4.
Reconstitution:	Centrifuge the vial before opening, reconstitute in sterile distilled water to a concentration of 0.1-1 mg/ml by gently pipetting 2-3 times, don't vortex.
Storage:	The lyophilized protein is stable at -20 °C for up to 1 year. After reconstitution, the protein solution is stable at -20 to -80 °C for 3 months or 1 week at 2 to 8 °C under sterile conditions. For extended storage, it is recommended to further dilute in working aliquots, avoid repeated freeze/thaw cycle.

Synonyms: CAF 1 antibody, CAF 1 subunit A antibody, CAF antibody, CAF I 150 kDa subunit antibody, CAF I p150 antibody, CAF-1 subunit A antibody, CAF-I 150 kDa subunit antibody, CAF-I p150 antibody, CAF1 antibody, CAF1A_HUMAN antibody, CAF1B antibody, CAF1P150 antibody, CHAF1A antibody

SDS-PAGE:



Safety Note: This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.