

Mouse/Rat CGRP (83-119) peptide

Catalogue No.:	P3543
Sequence:	SCNTATCVTHRLAGLLSRSGGVKDNFVPTNVGSEAF-NH ₂ (disulfide bridge Cys2-Cys7)
Molecular Formula:	C ₁₆₂ H ₂₆₂ N ₅₀ O ₅₂ S ₂
Formulation:	Lyophilized powder, TFA salt
Purity:	96.79%
Reconstitution:	Centrifuge the vial prior to opening, reconstitute in DMSO by gently pipetting 2-3 times, don't vortex
Storage:	The lyophilized powder is stable at -20°C for up to 1 year. Avoid repeated freeze/thaw cycle.
Safety note:	This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.