


Recombinant Human ARPC2

Catalogue No.:	P4099
Species:	Human
Uniprot ID:	O15144
Expression Region:	1-300
Host:	E.Coli
Tags:	N-terminal His Tag
Molecular Weight:	32.9 kDa under reducing conditions
Purity:	Greater than 90% as determined by SDS-PAGE.
Formulation:	Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 150 mM NaCl with 5% trehalose, pH7.4.
Reconstitution:	Centrifuge the vial before opening, reconstitute in sterile distilled water to a concentration of 0.1-1 mg/ml by gently pipetting 2-3 times, don't vortex.
Storage:	The lyophilized protein is stable at -20 °C for up to 1 year. After reconstitution, the protein solution is stable at -20 to -80 °C for 3 months or 1 week at 2 to 8 °C under sterile conditions. For extended storage, it is recommended to further dilute in working aliquots, avoid repeated freeze/thaw cycle.

Synonyms:

Actin related protein 2/3 complex subunit 2, Actin related protein 2/3 complex subunit 2 34kDa, Actin-related protein 2/3 complex subunit 2, ARC 34, ARC34, Arp2/3 complex 34 kDa subunit, ARP2/3 protein complex subunit 34, ARPC 2, Arpc2, ARPC2_HUMAN, p34 Arc, p34-ARC, p34Arc, PNAS 139, PNAS139, PRO 2446, PRO2446

SDS-PAGE:**Safety Note:**

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.