

Recombinant Human C4A

Catalogue No.:	P4644
Species:	Human
Uniprot ID:	P0C0L4
Expression Region:	1454-1744
Host:	E.Coli
Tags:	N-terminal His Tag or N-terminal His-IF2DI Tag, determined during production process
Molecular Weight:	TBD
Purity:	Greater than 90% as determined by SDS-PAGE
Formulation:	Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 150 mM NaCl with 5% trehalose, pH7.4
Reconstitution:	Centrifuge the vial prior to opening, reconstitute in sterile distilled water or buffer of interest to a concentration of 0.1-1mg/ml by gently pipetting 2-3 times, don't vortex.
Synonyms:	Acidic C4, Acidic complement C4, Basic C4, Basic complement C4, C3 and PZP-like alpha-2-macroglobulin domain-containing protein 2, C3 and PZP-like alpha-2-macroglobulin domain-containing protein 3, C4, Chido form, C4, Rodgers from, C4A, C4A anaphylatoxin, C4A2, C4A3, C4A4

- Storage:** The lyophilized protein is stable at -20°C for up to 1 year, the protein solution can be stored at 2-8°C for up to 1 week under sterile condition. For extended storage, it is recommended to further dilute in a buffer containing 0.1% BSA (make sure BSA will NOT interfere with your experiment) and store in working aliquots at -20°C to -80°C for 6 months. Avoid repeated freeze/thaw cycle. Avoid repeated freeze/thaw cycle.
- SDS-PAGE:** Test in progress
- Safety Note:** This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.