

Recombinant Human UPB1

Catalogue No.: P6900

Species: Human

Uniprot ID: Q9UBR1

Expression Region: 1-384

Host: E.Coli

Tags: N-terminal His Tag or N-terminal His-IF2DI Tag, determined during

production process

Molecular Weight: 42.1 kDa under reducing conditions

Purity: Greater than 90% as determined by SDS-PAGE.

Formulation: Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 150 mM NaCl

with 5% trehalose, pH7.4.

Reconstitution: Centrifuge the vial before opening, reconstitute in sterile distilled water to

a concentration of 0.1-1 mg/ml by gently pipetting 2-3 times, don't vortex.

Storage: The lyophilized protein is stable at -20 °C for up to 1 year. After

reconstitution, the protein solution is stable at -20 to -80 °C for 3 months

or 1 week at 2 to 8 °C under sterile conditions. For extended storage, it is

recommended to further dilute in working aliquots, avoid repeated


freeze/thaw cycle.



Synonyms:

AI195023, Beta alanine synthase, Beta ureidopropionase, BUP1, N carbamoyl beta alanine amidohydrolase, Ureidopropionase, Ureidopropionase, beta

SDS-PAGE:



Safety Note:

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.